

# University of Verona, School of Exercise and Sport Science, Laurea magistrale in Scienze motorie preventive ed adattate

Metodologia delle misure delle attività sportive

Thursday 24/10/2013

Luca P. Ardigò Ph.D.

## PA & ME

```
V'02, DLW issues
- ->ME;
- (mass?);
- specific activity efficiency? (PA/[ME-rME])
Luckily eff≈constant
- human species;
- walking main activity
HR issues
- ->ME (=kHR);
- <-external factors (e.g., stress, activity kind);
- latency
```

#### Observation issues

- simple, immediate;
- operator, video shooting;
- privacy?
- time taking (post-processing)

## Self-report diaries, questionnaires issues

- cheap;
- time taking (post-processing)
- correlation w/DLW higher than HR

## Common questionnaires

- Physical Activity Scale for the Elderly (PASE);
- Baecke questionnaire;
- Five-City questionnaire;
- Tecumseh questionnaire;
- Minnesota Leisure Time Physical Activity Questionnaire (MLTPA);
- Framingham questionnaire;
- Yale Physical Activity Survey (YPAS)

Ainslie et al., 2003

### Pedometry features

- waist;
- ->steps

#### Pedometer kinds

- electromechanical circuit based;
- electromagnetic circuit based;
- uniaxially accelerometric;
- ankle, shoe 1, 2 uni-, biaxially accelerometric


## Pedometry issues

- steps (i.e., most common PA kind m.u.) number;
- Japan standard Max e 3%


How many steps/day are enough?

- 10,000 (Hatano, 1993);
- Tudor-Locke et al., 2004:
- <5,000 sedentary lifestyle;
- 5,000 ÷ 7,499 typical daily activity that does not include exercise or sports and can be defined poorly active;
- 7,500 ÷ 9,999 includes a bit of extra-work (and/or fatiguing work) and can be defined a little active;
- -> 10,000 active lifestyle;
- -> 12,500 very active lifestyle

# Accuracy and precision



Good accuracy, poor trueness, poor precision



Low accuracy,
poor trueness,
good precision

# Accuracy and precision




Good accuracy, poor trueness, poor precision




Low accuracy,
poor trueness,
good precision

# Accuracy and precision



accuracy + precision = reliability

# Convergent validity




vs. accelerometer;

vs. observation;

vs. HR, V'O2, DLW;

vs. self-report diary

# Construct validity



w/age; w/anthropometry; w/fitness measures





FIGURE 1—Effect of speed on pedometer accuracy (percentage of actual steps) during treadmill walking.

stride #



(estimated) speed





FIGURE 1—Mean difference scores [(comparison – criterion pedometer)/criterion] ± SE as a percentage of the criterion estimated steps over a 24-h period. Positive difference scores represent overestimations, and negative difference scores indicate underestimations of steps compared with the criterion pedometer.

step/day #



400-m step #

# Pedometer accuracy/reliability

#### measures




FIGURE 2—Effect of BMI (25–29.9 kg·m $^{-2}$ , 30–35 kg·m $^{-2}$ , and >35 kg·m $^{-2}$ ) on the percent of actual steps recorded by the New-Lifestyles NL-2000 (NL) and Yamax Digiwalker SW-200 (SW). Error bars are standard deviation. \* Significantly different from actual steps; \* significantly different from the NL (P < 0.05).

(uniaxially accelerometric)

stride #

(electromechanical circuit based)

## Pedometer

## Final pedometry issues


- no discrimination of weight lifting, gradient legged locomotion, cycling, swimming, rowing;
- shoe or ankle accelerometric pedometer -> stride #

## Accelerometer

First generation

#### Rationale

$$a \leftarrow F (= m \times a) \leftarrow by paying ME$$





**Figure 5.** Metabolic energy expenditure (EE) is estimated by computing the signal magnitude area (SMA). The acceleration signal is rectified and then integrated. EE is then estimated by means of a linear regression.


Placement waist but also chest, back, wrist, ankle

## Accelerometer


- A piezo-electric or -resistive sensor;
- -> 'counts' number, intensity;
- no isometric force







## Accelerometer




Mathiee et al., 2004

gravitational acceleration acting along the sensitive axis.

Accelerometer

#### measures




**Figure 4.** Acceleration signal produced by a waist-mounted accelerometer aligned in the vertical (gravitational) direction, during a selection of basic daily movements. The acceleration signal is composed of the gravitational acceleration due to the posture of the subject and the acceleration due to body movement. g is the acceleration due to gravity, approximately 9.81 m s<sup>-2</sup>. The measured accelerations are dependent on the activity being performed. If the accelerometer was attached at a different point on the body, different acceleration signals would be recorded.

# Accelerometer

# 128 · 1). Raw -128 60 sec 128 2). Full-wave Rectified 60 sec 4000-3). Integrated (15 sec) 60 sec

FIGURE 2—Analytical processing of the acceleration data. 1. Raw: a 60-s window of a digitized raw signal collected at 32 Hz and using a 8-bit A/D conversion. 2. Rectification: all negative signal from (1) was turned into positive. 3. Integration: 15-s epochs.

## Accelerometer



## Accelerometer


```
- Uniaxial accelerometer -> -59% children DLW ME;
                          -59÷-50% old men DLW ME;
                          +50÷+60% old claudicants DLW ME;
                          -59% young women weekly DLW nME;
- Triaxial accelerometer -> +12÷+49% V'O2 locomotion ME;
                         -21÷-8% V'O2 gradient walking ME;
                         -68÷-53% V'O2 cycling ME;
                         -45÷-35% V'O2 daily activity ME;
                         -35% young women weekly DLW nME
```

# HR wrist monitor









#### HR measure issues

- HR <- environmental temperature and humidity, hydration status, posture, illness, stress, type of exercise (w/upper limbs or lower ones, continuous or intermittent), gender, age, body mass;
- w/≈3' latency

#### but...

- pedometer/accelerometer -> level legged locomotion;
   no exercise with upper limbs, walking and running on soft ground or on slopes, cycling, swimming, rowing;
- pedometer/uniaxial accelerometer -> no over 9 km/h running

#### HR measure issues

- (HR  $\geq$  90 bpm or  $\geq$  60% HRMax) ME = k HR;
- -30% daily ME;

## (partial) answer:

- FLEX HR method (Spurr et al., 1988): ME = k HR (subject, activity specific) use only @external load/HR > FLEX HR, i.e., average between maximum value during rest or sedentary activity, and minimum value during light activity;
- i.e., (HR < FLEX HR) ME = rME, (HR > FLEX HR) ME = k HR;
- -17÷+20% daily DLW ME



-> beat to beat recording -> HRV

# V'02 measure

- direct calorimetry in metabolic chamber;
- indirect calorimetry, respirometry @closed/open circuit -> V'O2, V'CO2, ...;
- < 8h





# V'02 measure








# V'02 measure



#### Option 1:

An operator can follow the swimmer by holding the K4b2 using a special rod (rod and harness are included in the standard packaging)



Option 2:

The K4b2 can be hung on a cable to be placed above the swimming pool lane